



QUANTUM[®]

#1 FOR REHAB POWER

Dedicated to Revolutionising the Power Chair Experience

TRU₄
BALANCE[®]
POWER POSITIONING SYSTEMS





WELCOME TO QUANTUM®

Quantum Rehab® was born out of the desire to delight customers with the most advanced, consumer-inspired complex rehab power wheelchairs and related technologies possible.

Clinically-focused and consumer-inspired, Quantum® products are engineered and manufactured to meet the widest range of complex rehab needs with functional, durable, stylised, high-performance designs to increase daily functionality and quality of life. At Quantum, consumer needs and wishes are the driving force. We're dedicated to not just meeting medical and clinical needs, but also quality-of-life needs. From the most advanced power seating for pressure management to USB ports, Bluetooth®, and fender lights, no consumer need is overlooked.

Our Australian headquarters is based in Melbourne with additional warehouses in New South Wales, Queensland and regional Victoria. All our Quantum powerchairs are built in our Melbourne office, servicing all of Australia.



TRU
BALANCE[®] 4
POWER POSITIONING SYSTEMS

TRU-Balance[®] 4 Power Positioning Systems provide a choice of 10 or 30 degrees of power anterior tilt.

Optional backup camera with 8.9 cm display can be mounted in multiple locations

Optional 10" power adjustable safe seat allows driving at speeds up to 5.6 km/h while elevated

STC (Smart Traction Control)



Q-Logic 3 Advanced Drive Controls with Bluetooth[®]

Stealth Products[®] cushion featuring exclusive Coolcore[®] technology

Front and rear LED fender lights for seeing and being seen (standard)



"RIDING IN THE 4FRONT 2 IS UNLIKE ANY CHAIR I HAVE EVER USED. THE 4FRONT 2 PUSHES THE LIMIT SET BY ITS PREDECESSOR BEYOND MY EXPECTATIONS."

CHRISTIAN BUDNEY

 PROUDLY BUILT IN AUSTRALIA

AVAILABLE WITH 10" SAFE SEAT ELEVATION



4FRONT[®] 2

FRONT-WHEEL DRIVE


FEATURING STC (SMART TRACTION CONTROL)

The 4Front[®] 2 is the next generation of the 4Front[®] Power Chair! New Smart Traction Control (STC), along with Smooth Ride Suspension (SRS), provide stability for smooth obstacle transitions and unmatched driving performance. Additional new features include 4-pole motors and redesigned caster arms for improved durability and aesthetics and reduced weight.

SPECS

- Maximum Speed: 10 k/h (5.6 km/h when seat is elevated)
- Weight Capacity: 136 kg
- Turning Radius: 65.8 cm
- Power Seat Elevation: Optional 25.4 cm (10") Safe Seat Elevation

COLOURS




Redesigned caster arms for improved durability

Automotive-grade SRS (Smooth Ride Suspension)

LED front and rear fender lights (front: white, rear: red)

 QUANTUM[®]
#1 FOR REHAB POWER



Fun Fact...

Stretto is the Italian word for *narrow*.



"TIGHT ROOMS, NARROW DOORWAYS, OUT-OF-REACH CABINETS AND VAN RAMP... NOW I CAN GO WHERE I WANT TO GO, ALL BECAUSE OF THE STRETTO!"

BRYSON FOSTER

OUR NARROWEST BASE IS ONLY
52 CM WIDE
WITH 12.5" DRIVE WHEELS

TRU BALANCE⁴
POWER POSITIONING SYSTEMS

TRU-Balance⁴ Power Positioning Systems provide a choice of 10 or 30 degrees of power anterior tilt.



Q-Logic 3 Advanced Drive Controls with Bluetooth[®] and USB charger

Stealth Products[®] cushion featuring exclusive Coolcore[®] technology

Front and rear LED fender lights

Optional 12" iLevel[®] power seat elevation allows driving at speeds up to 5.6 km/h while elevated

Toolless removal of rear shroud with battery access door

Independent SRS (Smooth Ride Suspension)

 PROUDLY BUILT IN AUSTRALIA

AVAILABLE WITH iLEVEL[®]
POWER SEAT ELEVATION

EDGE 3 STRETTO[®]

MID-WHEEL DRIVE

The Edge 3 Stretto[®] is the narrowest, most manoeuvrable complex rehab power base available in Australia! Aside from its narrow width, it is equipped with independent SRS (Smooth Ride Suspension), making it ideal for children, teenagers, and small adults. The Stretto is available with 31.8 cm (12.5") or 35.6 cm (14") drive wheels and optional 30.4 cm (12") of iLevel[®] power adjustable seat height.

SPECS

- Maximum Speed: 10 km/h (5.6 km/h when elevated with iLevel[®])
- Weight Capacity: 136 kg
- Turning Radius: 53.3 cm
- Power Seat Elevation: Optional 30.4 cm (12") iLevel

COLOURS



Q
QUANTUM[®]
#1 FOR REHAB POWER

Independent SRS (Smooth Ride Suspension)

The Stretto is the ideal pediatric power chair solution

Toolless access to electronics



TRU
BALANCE[®] 4
POWER POSITIONING SYSTEMS

TRU-Balance[®] 4 Power Positioning Systems provide a choice of 10 or 30 degrees of power anterior tilt

Optional 12" iLevel[®] power seat elevation allows driving at speeds up to 7.2 km/h while elevated

Featuring Extra Stability Technology[®] for safer transfers, reach and other activities of daily living



Q-Logic 3
Advanced Drive Controls with Bluetooth[®]

Stealth Products[®] cushion featuring exclusive Coolcore[®] technology

Upgraded SRS (Smooth Ride Suspension)



"THE EDGE 3 WITH ILEVEL HAS RESTORED MY INDEPENDENCE AT BOTH WORK AND HOME. ILEVEL HAS ALLOWED ME TO GET MORE INVOLVED AND START A DIALOGUE WITH PEOPLE, SOMETHING I WOULDN'T HAVE DONE BEFORE."

MERLISHA HENDERSON

 PROUDLY BUILT IN AUSTRALIA



EDGE[®] 3


MID-WHEEL DRIVE

The Edge[®] 3 offers upgraded SRS (Smooth Ride Suspension), which provides a smooth, comfortable ride and advanced stability. Optional iLevel[®] is available at 7.24 km/h with up to 30.4 cm (12") of power adjustable seat height.

SPECS

- Maximum Speed: 10 km/h (7.24 km/h when elevated with iLevel[®])
- Weight Capacity: 136 kg (See Q6 Edge[®] HD for increased capacity)
- Turning Radius: 52 cm
- Power Seat Elevation: Optional 30.4 cm (12") iLevel

COLOURS



USB charger (standard)

Front and rear LED fender lights for seeing and being seen (standard)

Drive wheel color accents (standard)

 **QUANTUM[®]**
#1 FOR REHAB POWER



TRUBALANCE³
POWER POSITIONING SYSTEMS

Q-Logic 3
Advanced Drive
Controls with
Bluetooth®

TRU-Balance® 3 Power
Positioning System

Optional 10" power seat
elevation allows for
driving at speeds up to
5.6 km/h while elevated

Featuring Extra Stability
Technology® for safer
transfers, reach and
other activities of daily
living



Stealth Products®
cushion featuring
exclusive Coolcore®
technology

Upgraded SRS
(Smooth Ride
Suspension)

NEWEST RELEASE

 PROUDLY BUILT IN
AUSTRALIA

AVAILABLE WITH

POWER SEAT ELEVATION

R-TRAK®


REAR-WHEEL DRIVE

The R-TRAK revolutionizes the rear-wheel drive power chair experience. Independent suspension on all four wheels combined with a compact base provide a smoother, more maneuverable driving experience. Standard features include fender lights, USB charger and drive wheel color accents. The R-TRAK is also available with our 25.4 cm (10") safe seat elevation, allowing the client to drive up to 5.6 km/h while elevated.

SPECS

- Maximum Speed: 10 km/h
- Weight Capacity: 136 kg
- Turning Radius: 65 cm
- Power Seat Elevation: Optional 25.4 cm (10") elevation

COLOURS



Advanced Drive Controls
with Bluetooth®

Superior, automotive-grade SRS
(Smooth Ride Suspension)

Unmatched driving
performance and ride quality


QUANTUM®
#1 FOR REHAB POWER





Q6 EDGE® 2.0

MID-WHEEL DRIVE

SPECS

- Maximum Speed: 10 km/h (5.6 km/h when elevated)
- Weight Capacity: 136 kg
- Turning Radius: 52 cm (with 15.2 cm castors)
- Power Seat Elevation: Optional 30.4 cm iLevel®

COLOURS



4 FRONT 2 HD

MID-WHEEL DRIVE

SPECS

- Maximum Speed: 8 km/h
- Weight Capacity: 204 kg
- Turning Radius: 65.7 cm

COLOURS



Q6 EDGE® Z

MID-WHEEL DRIVE

SPECS

- Maximum Speed: 10 km/h (5.6 km/h when elevated)
- Weight Capacity: 160kg
- Turning Radius: 50 cm
- Power Seat Elevation: Optional 30.4 cm iLevel®

COLOURS



Q6 EDGE® HD

BARIATRIC MID-WHEEL DRIVE

SPECS

- Maximum Speed: 7.5 km/h (5.6 km/h when elevated)
- Weight Capacity: 204.1 kg (181.4 kg with iLevel®)
- Turning Radius: 50 cm
- Power Seat Elevation: Optional 25.4 cm iLevel®

COLOURS




Q1450

BARIATRIC FRONT-WHEEL DRIVE

SPECS

- Maximum Speed: 7.2 km/h
- Weight Capacity: 272 kg
- Turning Radius: 67.3 cm

COLOURS




J623 2.0

MID-WHEEL DRIVE

SPECS

- Maximum Speed: 10 km/h
- Weight Capacity: 136 kg
- Turning Radius: 52.1 cm

COLOURS



* ADDITIONAL J623 2.0 COLOURS INCLUDE: CARBON FIBRE, EMERALD GREEN, ORANGE CRUSH, SILVER BULLET, & YELLOW BLAZE



QUANTUM
iLevel
power chairs

TRU-Balance® 4 Power Positioning Systems

WITH ANTERIOR TILT AND MEMORY SEATING

Maximum weight capacity of up to 136 kg

Up to 8 customizable memory seating positions recall the saved positions of all actuators, offering the choice between sequential or synced movement for each seat function.

Seat elevation combined with anterior tilt provides optimal access.

Consumers may find it easier to perform ADLs such as: accessing sinks and cooktop surfaces, meal preparation and reaching items in cabinets.



TRU-Balance® 4 with Dual Actuator AFP shown on the Edge 3 Stretto®

Chest Bar

DUAL ACTUATOR POWER ARTICULATING FOOT PLATFORM

Lowers to the ground to assist with transfers

Articulation range of 0-7.25"

Elevation range of 85-165°

Only available on TRU-Balance® 4



Knee Blocks

TRU-BALANCE® 4

POWER POSITIONING SYSTEMS

FEATURES

- Anterior tilt options: 10 and 30 degrees
- Up to eight customizable memory seating positions
- Chest bar and knee blocks available (knee blocks and chest support required with 20 or 30 degrees of anterior tilt)
- Available on the Edge 3 Stretto®, 4Front® 2, and Edge® 3
- Optional Accu-Trac Advanced Tracking Technology
- Optional WC19 crash-test approved occupied transit kit option
- Available seat sizes: 16-22" width and 14-22" depth

POWER RECLINE

85° - 165° of recline

POWER POSTERIOR TILT

0° - 50° of tilt

Maintains center of gravity shift

DRIVE SPEEDS IN ANTERIOR TILT

Up to 10° of anterior tilt, up to full lift: 3.2 km/h

10-20° of anterior tilt, up to 15.3 cm (6") of lift: 0.8 km/h

10-20° of anterior tilt, over 15.3 cm (6") of lift: drive lockout

Over 20° of anterior tilt: drive lockout

TRU-Balance® 4 Power Recline



TRU-Balance® 4 Power Posterior Tilt



TRU-Balance® 4 with 12" iLevel® Power Elevation





"COOKING IS SO MUCH EASIER NOW THAT I CAN REACH THINGS. I CAN'T IMAGINE MY LIFE WITHOUT ILEVEL."

STEPHANIE WOODWARD

I LEVEL®

POWER SEAT ELEVATION

Standing Height + Speed + Stability = Authentic iLevel (like no other)

iLevel® power seat elevation technology features a revolutionary new power chair design for a better experience. iLevel improves daily functionality and quality of life by allowing operation of the power chair with the seat elevated while at walking speed.

FEATURES

- 30.4 cm (12") of power adjustable seat height while at walking speed; maximum of 25.4 cm (10") height on the Q6 Edge® HD
- iLevel's patented Extra Stability Technology® enhances safety for transfers, reach, and numerous other activities of daily living
- Available on the Edge® 3 and Q6 Edge® 2.0
- Front and rear LED fender lights come standard with iLevel®



Safer Transfers



Improved Reach



Social interaction at eye level



TRU-Balance® 3 Power Adjustable Seat Lift with iLevel® technology

Maximum weight capacity of up to 136 kg

Provides extra stability to ensure maximum safety for transfers, reach, and numerous other activities of daily living

Maximum speed of 7.24 km/h while elevated at 30.4 cm of lift on the Edge 3

Seat-to-floor-height as low as 44.45 cm on the Edge 3

LED fender lights included



POWER ARTICULATING FOOT PLATFORM

Broad articulation range for enhanced adjustment

204 kg weight capacity

Unmatched adjustment, durability, and strength

Multiple calf pads, calf panel, and footplate options

TRU-BALANCE® 3

POWER POSITIONING SYSTEMS

FEATURES

- Fine adjustments can be made even while client is in the chair
- Seat width can be ordered independently of the back width at no additional charge
- Power tilt, power recline, and iLevel® are retrofittable

POWER RECLINE

85° - 175° of recline Available in matching back shroud colours

136 kg weight capacity

Available with power tilt and power adjustable seat lift

POWER TILT

50° of tilt Maintains center of gravity shift

Lower seat-to-floor height

Unmatched adjustment and durability

HEAVY DUTY POWER TILT & RECLINE

204 kg weight capacity Available in matching back shroud colours

Up to 50° of tilt 85° - 175° of recline

TRU-Balance® 3 Power Recline



TRU-Balance® 3 Power Tilt



TRU-Balance® 3 Heavy Duty Power Tilt and Recline





MORGAN STEWARD ENJOYS THE SPEED, ELEVATION, AND RECLINE OF HIS ILEVEL® POWER CHAIR.

MORGAN STEWARD

Quantum® paediatric seating with iLevel® 30.4 cm power seat elevation

SEATING

AND POSITIONING

SYNERGY® SEATING

- Easily adjustable, this versatile design can grow with the user
- Several width/depth adjustment packages ranging from 10" to 32"
- Accepts a wide variety of after-market seats and backs
- Heavy duty option provides a weight capacity of 272 kg and seating up to 32" x 28" (model specific)



Q6 Edge® 2.0 with Synergy® seating

BARIATRIC TILT

- TRU-Balance HD Power Positioning Systems offers increased weight capacities and expanded seating sizes
- Configurations with a 204kg weight capacity include HD power tilt and HD power recline
- 10" power adjustable seat lift is available with a 181kg weight capacity on the Q6 Edge® HD the 4Front® 2 HD.
- Bariatric Tilt for configurations up to 272kgs is also available on the Q1450.

Synergy® Seating



Contoured Seating



Stealth Products® CUSHIONS

FEATURING EXCLUSIVE
coolcore
CHEMICAL-FREE COOLING TECHNOLOGY



30+ UPF



Dry




Chemical Free



Breathable Comfort

"While the other makers of cooling fabrics use a number of chemicals to enhance wicking, transportation and evaporation, Coolcore fabrics are chemical free. To make that possible, we changed the essence of cooling fabric technology right down to its very core."

COOLCORE®



STEALTH PRODUCTS®

CUSHIONS

Stealth Products® Cushions feature patented Coolcore® technology, a chemical free, four-way stretch material that offers moisture wicking, moisture transportation, and moisture evaporation properties. Our cushions also use a colour guide that identifies the suggested function family by a unique patch.

FEATURES

- Chemical free
- Regulates temperature
- Manages and distributes your body heat



Convenient built-in side/front pouches









Constructed with FRG Sunmate flame-retardant material

Easy-pull rubber loop zipper handles

Q
QUANTUM®
#1 FOR REHAB POWER





	TRU-Comfort 2 SPP	Essence SPP	Solution® SPP	Spectrum® Foam/Gel SPP	Premiere P	Glacial SP	Zen SP	Simplicity G
DESCRIPTION¹	 <p>Multi-layered foam Deep contour, CoolCore® cushion cover has "pockets" on sides</p>	 <p>Low profile Tru-Comfort 2, lateral and pre-ischial shelf Contoured, HD molded foam base, coccyx and sacral relief</p>	 <p>Standard viscoelastic foam top with twin cell aqueous viscous gel (heavier gel) insert in ischial well, mid-level contour</p>	 <p>SunMate FRG foam top (SunMate FRG XXSoft (blue); Medsoft (green)) with multi-layered foam insert or twin cell ultra-lightweight viscous gel insert in ischial well, sacral relief</p>	 <p>Gray X-Soft SunMate FRG foam topper, high density positioning wedges built in for offloading.</p>	 <p>New SIM G base, fluid resistant overlay, top and sides TPU heat-sealed, low profile, lightweight cooling gel topper</p>	 <p>New SIM G base, SunMate FRG XXSoft (blue topper), low profile, lightweight</p>	 <p>High-density foam Stealth-Tek cover standard More contour to sacral/coccyx relief well</p>
WEIGHT CAPACITY	14-15" W: 79.4 kg 16-21" W: 158.8 kg 22+" W: 204.1 kg	10-15" W: 79.4 kg 16-21" W: 158.8 kg 22" W: 204.1 kg	10-15" W: 79.4 kg 16-20" W: 158.8 kg 22+" W: 226.8 kg	14" W: 79.4 kg Foam & Gel 16-20" W: 158.8 kg Foam & Gel 22+" W: 294.8 kg Foam; 550 lbs. Gel	14-15" W: 79.4 kg 16-21" W: 158.8 kg 22" W: 226.8 kg	10-15" W: 79.4 kg 16-21" W: 158.8 kg 22+" W: 181.4 kg	10-15" W: 79.4 kg 16-21" W: 158.8 kg 22+" W: 181.4 kg	10-15" W: 79.4 kg 16-21" W: 158.8 kg 22+" W: 181.4 kg
OVERALL WEIGHT	1.8 kg	1.6 - 2.3 kg	2.4 kg	Foam: 1.4 kg Gel: 1.6 kg	1.6 kg	1.6 kg	1.1 kg	.9 kg
HEIGHT	11.4 cm / 4.5"	8.9 cm / 3.5"	7.6 cm / 3"	7.6 cm / 3"	8.9 cm / 3.5"	6.4 cm / 2.5"	7.6 cm / 3"	5 cm / 2"
DIMENSIONS²	14x14" up to 24x24"	10x10" up to 22x24"	10x10" up to 24x24"	Foam: 14x18" up to 32x28" Gel: 14x18" up to 24x26"	14x13" up to 22x22"	10x10" up to 24x24" Odd sizes in 15-21" W	10x10" up to 24x24" Odd sizes in 15-21" W	10x10" up to 24x24" Odd sizes in 15-21" W

1) CoolCore® is standard on all cushions except the Simplicity G.

ADI Q-Back 2.0



The ADI Q-Back 2.0 from Stealth Products® is available on TRU-Balance® 3 Seating Systems and is offered in three configurations:

- Active Contour ADI Q-Back 2.0
- Deep Contour ADI Q-Back 2.0
- Extra Deep Contour ADI Q-Back 2.0
- Coolcore® cover comes standard

Sizes for static/power tilt/power recline seating:

- Support heights: 13", 16", 18", 20", 22"
- Widths: 14", 15", 16", 17", 18", 19", 20", 22"
- Additional back sizes are available with special order

3" OF BACK SUPPORT HEIGHT ADJUSTMENT BUILT INTO THE TRU-BALANCE® 3 BACK SPINE.

TRU-Comfort 2 Back




- Adjustable positioning and pressure management
- Compatible with both TRU-Balance®3 Positioning Systems and the frame system
- Four-way stretch, recovery, fluid-resistant, breathable cover

UniLink™ Hardware System

- New system from Stealth Products® for mounting thigh guides and elbow blocks onto the TRU-Balance® 3 Seating System
- Elbow blocks in 3"x10", 3"x3", and 4"x4" pad sizes are made with patented Coolcore® technology cover material
- Thigh guides are available in 4"x4", 4"x6", 4"x8", 4"x10", and 4"x12" pad sizes

THIGH GUIDE OPTIONS





"MY EDGE GIVES ME THE OPPORTUNITY TO USE ADVANCED TECHNOLOGY THAT GIVES ME SO MUCH INDEPENDENCE."

ZOE HERNANDEZ



Q-LOGIC 3

ADVANCED DRIVE CONTROLS


The Q-Logic 3 Advanced Drive Controls System from Quantum Rehab® is like no other power wheelchair electronics system. Designed from clinician and consumer feedback, Q-Logic 3 is user-friendly and highly customisable.

FEATURES


- Wireless Bluetooth® programming through Windows devices
- Bluetooth access to computer, tablet, and phone functions
- iAccess programmable push button and toggle options for customised seating access

Clinic Mode:

- Quickly set up multiple input devices in client evaluations
- Automatically assigned drive modes to input devices



Q-Logic 3 home screen



Q-Logic 3 drive screen

ENHANCED DISPLAY

The optional Q-Logic 3 EX Enhanced Display serves as a separate screen used with specialty drive controls to display system data. With the Enhanced Display's built-in environmental controls, Q-Logic 3 can operate a wide range of appliances that use infrared (IR) remote controls, such as televisions, and serves as a Bluetooth® computer mouse emulator.



Q-LOGIC 3 EX ENHANCED DISPLAY

Q-LOGIC 3e JOYSTICK CONTROLLER

- Industry-first, magnetic charging port with adapter can be used with any current XLR charger
- Magnetic wall holder keeps charger plug within reach and off floor
- Integrated puddle light provides greater visibility



Q-LOGIC 3 EX JOYSTICK CONTROLLER

Q-Logic 3 iAccess Programmable Switchbox



Q-Logic 3 Stand Alone Joystick



Q-Logic 3 Specialty Controls Interface Module





ACCESSORIES

VR2 6 KEY JOYSTICK



JOYSTICK OPTIONS



SIP-N-PUFF



Stealth® i-Drive® Sip-N-Puff

MULTI-SWITCH



Stealth® i-Drive® 4 Stealth Pro, 3-Switch Head Array



Stealth® i-Drive® Tri-Array, 3-Switch Head Array



Stealth® i-Drive® Curved Array

Q-LOGIC 3



Q-Logic 3 EX Joystick



Q-Logic 3e Joystick



Q-Logic 3 EX Enhanced Display



Q-Logic 3 Stand Alone Joystick



Q-Logic 3 iAccess Programmable Switchbox



Q-Logic 3 Specialty Interface Module

SINGLE SWITCH



Stealth® Egg Switch



Mo-Vis™ Twister Basic Button



Mo-Vis™ Twister Satellite Button

TRAY ARRAYS



Stealth® i-Drive® Eclipse Tray with Four Capacitive Proximity Sensors



Stealth® i-Drive® Eclipse Tray with Four Fiber Optic Proximity Sensors

PROPORTIONAL DRIVE CONTROLS



Mo-Vis™ All-Round Joystick



Mo-Vis™ All-Round Light Joystick



Mo-Vis™ Multi Joystick



Stealth® i-Drive® Control Harness



Stealth® Precision Mini Proportional Joystick



Stealth® Mushroom Joystick



Mo-Vis™ Micro Joystick

MISC. ACCESSORIES



Front and Rear LED Fender Lights



Hydration System



Transfer Bars (Pair)



Quantum® Backpack



Personal Item Hook



Phone Holder





















Reversing Camera



Glove Box

COLOUR GUIDE

		STRETTO®	EDGE 3	Q6 EDGE® 2.0	4FRONT® 2	Q6 EDGE® Z	Q6 EDGE® HD	Q1450	J623 2.0*	R-TRAK®
	Deep Purple	•	•	•	•	•	•			•
	Sugar Plum								•	
	Midnight Blue		•	•	•	•	•			•
	Blue Jean	•	•							
	Electric Blue			•	•	•	•	•	•	•
	Military Green		•							
	Night Train		•	•	•					
	Back in Black	•	•	•	•	•	•	•	•	•
	Grey Street			•						
	White Satin			•						
	Green Machine	•	•	•		•	•			•
	Lemon Crush		•	•						•
	Evolution Orange		•	•	•					•
	Red Corvette	•	•	•	•	•	•	•	•	•
	Cherry Bomb		•							
	Raspberry Sorbet	•	•							
	Pretty in Pink			•					•	•
	Lime Light				•				•	

* ADDITIONAL J623 2.0 COLOURS INCLUDE: CARBON FIBRE, EMERALD GREEN, ORANGE CRUSH, SILVER BULLET, & YELLOW BLAZE. CONTACT YOUR NEAREST AUTHORISED QUANTUM PROVIDER TO LEARN MORE.





	EDGE 3 STRETTO®	EDGE® 3	Q6 EDGE® 2.0
DRIVE CONFIGURATION	Mid-Wheel Drive	Mid-Wheel Drive	Mid-Wheel Drive
DRIVE WHEELS	35.5 cm (14") or 31.7 cm (12.5") Pneumatic or Flat Free Tyres	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²
CASTOR WHEELS	Front: 12.7 cm (5") or 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid	Front: 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid	Front: 15.2 cm (6") Solid ¹ 20.3 cm (8") Solid ² Rear: 15.2 cm (6") Solid ¹ 20.3 cm (8") Solid ²
ANTI-TIP WHEELS	N/A	N/A	N/A
MAXIMUM SPEED³	Up to 10 km/h (5.6 km/h w/ iLevel) ¹⁰	Up to 10 km/h (7.2 km/h with iLevel)	Up to 10 km/h (5.6 km/h w/ iLevel) ¹⁰
GROUND CLEARANCE	7.3 cm (2.9")	6.8 cm (2.7") (battery tray)	6.8 cm (2.7") (battery tray)
TURNING RADIUS⁴	53.3 cm (21")	52 cm (20.5")	52 cm (20.5") w/ 15.2 cm (6") castors 59 cm (23.2") w/ 20.3 cm (8") castors
OVERALL LENGTH⁸	82.5 cm (32.5")	81.8 cm (32.2") without foot riggings	90 cm (35.5") w/ 15.2 cm (6") castors 102.6 cm (40.4") w/ 20.3 cm (8") castors
BASE WIDTH	55.2 cm (21.75") (35.5 cm (14") drive wheels) 51.9 cm (20.47") (31.5 cm (12.5") drive wheels)	61.9 cm (24.4")	61.9 cm (24.4")
SEAT SIZES	TRU-Balance® 3 and 4	TRU-Balance® 3 and 4	TRU-Balance® 3
TRU-BALANCE® 4	W: 16-22" D: 14-22"	W: 16-22" D: 14-22"	N/A
TRU-BALANCE 3® STATIC	W: 12-20" D: 12-20"	W: 12-22" D: 10-22"	W: 12-22" D: 12-22"
TRU-BALANCE® 3 POWER TILT	W: 12-20" D: 12-20"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"
TRU-BALANCE® 3 LIFT & TILT	W: 12-20" D: 12-20"	W: 12-22" D: 12-22"	W: 12-22" D: 12-22"
SEAT-TO-FLOOR HEIGHTS^{13,14}	TRU-Balance® 3 and 4	TRU-Balance® 3 and 4	TRU-Balance® 3
TRU-BALANCE® 4	17.5-19.5"	17.5-19.5"	17.5-19.5"
SYNERGY / STATIC	16-18"	16.5-18.5"	16.5-18.5"
TRU-BALANCE® 3 POWER TILT⁷	17.5-19.5"	17.5-19.5"	17.5-19.5"
TRU-BALANCE® 3 LIFT & TILT⁷	17.5-19.5"	17.5-19.5"	17.5-19.5"
BATTERY SIZE⁵	MK 48Ah(2)	MK 50Ah(2)	MK 50Ah(2)
BATTERY WEIGHT¹¹	17 kg each	17.2 kg each	17.2 kg each
AVAILABLE ELECTRONICS⁹	120A Q-Logic 3	120A Q-Logic 3	120A Q-Logic 3
BATTERY CHARGER	8A Off-Board	8A Off-Board	8A Off-Board
MOTOR PACKAGES	2-Pole	2-Pole or 4-Pole	2-Pole
WEIGHT CAPACITY⁶	136 kg	136 kg	136 kg
BASE WEIGHT	57.3 kg	62.8 kg	62.8 kg w/ 15.2 cm (6") castors 67.1 kg w/ 20.3 cm (8") castors

	4 FRONT® 2	4 FRONT® 2 HD	Q6 EDGE® Z	Q6 EDGE® HD ¹⁵
DRIVE CONFIGURATION	Front-Wheel Drive	Front-Wheel Drive	Mid-Wheel Drive	Mid-Wheel Drive
DRIVE WHEELS	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flate Free Tyres ²
CASTOR WHEELS	Rear: 22.8 cm (9") Solid	Rear: 22.8 cm (9") Solid	Front: 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid	Front: 15.2 cm (6") Solid Rear: 15.2 cm (6") Solid
ANTI-TIP WHEELS	Front: 10 cm (4") Solid	Front: 10.1 cm (4") S	N/A	N/A
MAXIMUM SPEED³	Up to 10 km/h (5.6 km/h w/ elevation)	Up to 8 km/h (4.3 km/h w/ elevation)	Up to 10 km/h (5.6 km/h w/ iLevel)	Up to 7.6 km/h (5.6 km/h w/ iLevel)
GROUND CLEARANCE	7.6 cm (3") (battery tray)	7.6 cm (3") (battery tray)	7.6 cm (3") (frame base)	7.6 cm (3") (frame centre)
TURNING RADIUS⁴	65.7 cm (25.9")	65.7 cm (25.9")	50 cm (19.7")	50 cm (19.7")
OVERALL LENGTH⁸	101 cm (39.9")	101 cm (39.9")	90.4 cm (35.6")	90.4 cm (35.6")
BASE WIDTH	60.9 cm (24")	60.9 cm (24")	66.5 cm (26.2")	66.5 cm (26.2")
SEAT SIZES	TRU-Balance®3 and 4	TRU-Balance®3	TRU-Balance® 3	TB2 HD TB3 HD
TRU-BALANCE® 4	W: 16-22" D: 14-22"	N/A	N/A	N/A N/A
TRU-BALANCE 3® STATIC	W: 12-22" D: 12-22"	HD W: 20-24" D: 20-24"	W: 12-24" D: 12-24"	W: 16-28" D: 20-28" W: 16-32" D: 20-28"
TRU-BALANCE® 3 POWER TILT	W: 12-22" D: 12-22"	HD W: 20-24" D: 20-24"	W: 12-24" D: 14-24"	W: 16-28" D: 20-28" W: 16-28" D: 20-25"
TRU-BALANCE® 3 LIFT & TILT	W: 12-22" D: 12-22"	HD W: 20-24" D: 20-24"	W: 12-24" D: 10-24"	N/A W: 16-24" D: 20-25"
SEAT-TO-FLOOR HEIGHTS^{13,14}	TRU-Balance® 3 and 4	TRU-Balance®3	TRU-Balance® 3	TB2 HD TB3 HD
TRU-BALANCE® 4	17.75-19.75"	N/A	18-20"	N/A 17.7-19.7"
SYNERGY / STATIC	15.2-17.2"	15-17"	16.5-18.5"	16.5-18.5" 17.5-19.5"
TRU-BALANCE® 3 POWER TILT⁷	17.75-19.75"	17.75-19.75"	17.5-19.5"	18.2-20.9" 18.5-20.5"
TRU-BALANCE® 3 LIFT & TILT⁷	17.75-19.75"	17.75-19.75"	18-20"	N/A 17.7-19.7"
BATTERY SIZE⁵	MK 70Ah(2)	MK 70Ah(2)	MK 70Ah(2)	MK 70Ah(2)
BATTERY WEIGHT¹¹	23.3 kg each	23.3 kg each	23.5 kg each	23.5 kg each
AVAILABLE ELECTRONICS⁹	120A Q-Logic 3	120A Q-Logic 3	120A Q-Logic 3	120A Q-Logic 3
BATTERY CHARGER	8A Off-Board	8A Off-Board	8A Off-Board	8A Off-Board
MOTOR PACKAGES	4-Pole	4-Pole	4-Pole	4-Pole ¹
WEIGHT CAPACITY⁶	136 kg	204 kg	160 kg	204.1 kg ¹²
BASE WEIGHT	58.9 k	58.9 kg	65 kg	65 kg

1 Standard. 2 Optional. 3 Speed varies with user weight, terrain type, battery charge, battery condition and tyre pressure. 4 Without front riggings. 5 Battery capacity may change when adding power positioning systems. 6 Weight capacity is dictated by the selection of seating. 7 All measurements +/- 1/8" with 18" x18" seat, fixed towers, 136 kg weight capacity and solid tyres (add 3/8" with pneumatic tyres). Depths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in their lowest position; does not include elevating system unless specified.

8 Without foot riggings. 9 All Quantum units can be ordered with speciality controls. 10 10 km/h speed availability only with Q-Logic 3. 11 Battery weight may vary +/- 1.3 kg based on manufacturer. 12 Weight capacity is reduced to 181 kg with iLevel® power adjustable seat height. 13 Seat-to-floor heights are measured with the unit unloaded. Due to the 4Front's suspension system, seat-to-floor heights are approximately 0.6 cm lower when the unit is loaded. 14 41.9 cm seat-to-floor height is not available with power tilt or iLevel®. 49.5 cm seat-to-floor height is not available with static seat



	QUANTUM® 1450	J623 2.0	R-TRAK
DRIVE CONFIGURATION	Front-Wheel Drive	Mid-Wheel Drive	Rear-Wheel Drive
DRIVE WHEELS	35.5 x 10 cm (14x4" wide) Flat Free Tyres	35.5 cm (14") Pneumatic or Flat Free Tyres 14x4" wide Flat Free Tyres ²	35.6 cm (14") Pneumatic or Solid Tyres 14x4" wide Flat Free Tyres ²
CASTOR WHEELS	22.8 cm (9") Solid	Front: 15.2 cm (6") Solid ¹ 20.3 cm (8") Solid ² Rear: 15.2 cm (6") Solid ¹ 20.3 cm (8") Solid ²	Front: 22.9 cm (9") Solid or Pneumatic
ANTI-TIP WHEELS	Front: 7.5 cm (3") Solid	N/A	Rear: 10 cm (4") Solid
MAXIMUM SPEED³	Up to 7.2 km/h	Up to 10 km/h (2.5 km/h w/ elevation)	Up to 10 km/h (5.6 km/h w/ elevation)
GROUND CLEARANCE	5 cm (2") (anti-tips)	6.9 cm (2.7") (battery tray)	10 cm (3.97") (battery tray)
TURNING RADIUS⁴	67.3 cm (26.5")	52 cm (20.5")	65.3 cm (25.7")
OVERALL LENGTH⁸	90.4 cm (35.6")	90.2 cm (35.5") without foot riggings	88.3 cm (34.8")
BASE WIDTH	74 cm (29.2")	61.9 cm (24.4")	61.4 cm (24.2")
SEAT SIZES	TRU-Balance® 2 HD Tilt	TRU-Balance® 3	TRU-Balance® 3
TRU-BALANCE® 4	N/A	N/A	N/A
TRU-BALANCE 3® STATIC	W: 20-32" D: 20-28"	W: 12- 22" D: 12-22"	W: 12- 22" D: 12-22"
TRU-BALANCE® 3 POWER TILT	W: 20-32" D: 20-28"	W: 12- 22" D: 12-22"	W: 12- 22" D: 12-22"
TRU-BALANCE® 3 LIFT & TILT	N/A	N/A	W: 12- 22" D: 12-22"
SEAT-TO-FLOOR HEIGHTS^{13,14}	TRU-Balance® 2 HD Tilt	TRU-Balance® 3	TRU-Balance® 3
TRU-BALANCE® 4	N/A	17.5-19.5"	17.2-19.2"
SYNERGY / STATIC	16.5-18.5"	16.5-18.5"	16-18"
TRU-BALANCE® 3 POWER TILT⁷	18-20"	17.5-19.5"	17.2-19.2"
TRU-BALANCE® 3 LIFT & TILT⁷	N/A	17.5-19.5"	17.2-19.2"
BATTERY SIZE⁵	MK 70Ah(2)	MK 50Ah(2)	MK 70Ah(2)
BATTERY WEIGHT¹¹	23.5 kg each	17.2 kg each	23.3 kg each
AVAILABLE ELECTRONICS⁹	90A VR2 6-Key 120A Q-Logic 3	90A VR2 6-Key	120A Q-Logic 3
BATTERY CHARGER	8A Off-Board	8A Off-Board	8A Off-Board
MOTOR PACKAGES	4-Pole	2-Pole	4-Pole
WEIGHT CAPACITY⁶	272 kg	136 kg	136 kg
BASE WEIGHT	88.4 kg	62.8 kg	58 kg



or power recline. 15 The Q6 Edge® HD is not available in Lime Light NOTE: The information contained herein is correct at the time of publication.



SERVICE AND SUPPORT



QUANTUM®
AUSTRALIA


**WEBSITE
AND
PORTAL**

**COMING
SOON**



SOCIAL MEDIA

Stay up to date on the latest Quantum® events and product releases through Twitter, Facebook, YouTube, and Instagram.




QUANTUMREHAB.COM

Quantum's official website is the ultimate source for everything you need to know about Quantum® products. Highlights include detailed product specifications, order forms, videos, user testimonials, and much more.


LIFE.QUANTUMREHAB.COM

Your one-stop resource for articles and info on Quantum's® life-changing power wheelchairs and complex rehab technologies.



Q3DC.QUANTUMREHAB.COM (ECON-W)

Quantum Rehab's® login page for the Q-Logic 3 cloud-based wireless programming. Q-Logic 3 is our newest generation of electronics available with improved features and benefits inspired by consumer, provider and therapist feedback.



INTERACTIVE ASSIST

Q-Logic 3 now features Interactive Assist. Interactive Assist from Quantum Rehab® is an app which provides a direct real-time connection from a power chair's electronics to a provider's programming station. The connection provides the technician with full system and diagnostics information along with tools, such as a real-time mirror image of the power chair's electronics display, to eliminate wasteful trips.

- Change Log History
- Saving XPA Files
- Motor Usage History
- Text Chatting
- Programming IR Menu
- Access to Bluetooth Features
- Dealer-to-Dealer Programming Access
- Screen Mirroring
- Real-time Monitoring
- Efficiency Improvements
- Troubleshooting
- Battery Management
- Diagnostic History
- Basic ProgramminG
- Remote Access





Download Interactive Assist today!

ECON - I

The Econ-I app is now available for download on Apple® iPhone® devices, allowing you to program Quantum® Power Chairs conveniently from your iPhone or iPad. Compatible devices include the iPhone 6S/6s Plus/SE or higher, and the iPad Air/iPad Mini 2/iPad Pro/iPad 5th generation or higher. All devices must support iOS version 12 or higher.



Download ECON - I today!




PQ365 APP

ALL QUANTUM® PRODUCT INFORMATION IN ONE APP!

Designed for providers, the free, easy-to-use app helps busy professionals stay up to date with the latest information and allows access all relevant content in one quick and convenient place.

The app has many cool features, including, serial and order number search, order delivery status, invoice number look-up, a voice search function, an image share option, technical support, marketing material and educational resources. Additional features include a quick-search option, 360-degree product views, push notifications, news, events, specifications sheets, brochures, owner's manuals, basic operating instructions and more.



DOWNLOAD PQ-365 TODAY!



STEALTH RESOURCES

STEALTHPRODUCTS.COM



Stealth Products® offers a range of versatile products, including head and body support/positioning, joystick mounts, trays, seating, cushions, and paediatric mobility products.



CUSHIONS.STEALTHPRODUCTS.COM/

Learn about our Stealth Products range of cushions featuring our exclusive Coolcore® Technology.



Scan here to access the Stealth Cushions Catalogue



ADI.STEALTHPRODUCTS.COM/

Learn about our Stealth Products range of ergonomically designed ADI® back supports and seating systems.



Scan here to access the ADI Catalogue



UNILINK.COM

Learn about our Stealth Products UniLink range of strong, unique, and simplified components.



Scan here to access the UniLink Catalogue



HEADSOLUTIONS.COM

Learn about our Stealth Products Head Solutions range of versatile, high quality head and neck positioning systems and mount options.



Scan here to access the Head Solutions Catalogue



QUANTUM®

1 FOR REHAB POWER



Scan here to access
our Australian website