

Q-LOGIC 2

DRIVE CONTROL SYSTEM



Connecting Lives

 @QuantumRehab

 youtube.com/QuantumRehab

 facebook.com/QuantumRehab



20-24 Apollo Drive, Hallam, Victoria 3803
www.pridemobility.com.au • www.quantumrehab.com
Tel: 03 8770 9600 • Fax: 03 9703 2720



AU-QLGC-11-10-14

www.quantumrehab.com

Q-LOGIC 2

DRIVE CONTROL SYSTEM

The Q-Logic 2 Drive Control System by Quantum Rehab® is the most advanced power chair electronics package available. The Q-Logic Series creates an adaptable, expandable control system that serves even the most complex needs while delivering a host of exclusive features.

FEATURES

- More programming and customization options for a smooth ride
- Larger 2.8" LCD screen
- Built-in ambient light sensor automatically adjusts screen brightness
- Controller armor available in color matching options
- Enhanced horn volume
- Bluetooth access to computer, tablet and phone functions

Q-Logic 2 EX Enhanced Display



The optional Enhanced Display serves as a separate screen used with specialty drive controls to display system data. With the Enhanced Display's built-in environmental controls, Q-Logic 2 can operate a wide range of appliances that use infrared (IR) remote controls, such as televisions, and serves as a Bluetooth computer mouse emulator.

Optional Modules



Programming Devices



Q-Logic 2 Joystick

Non-Expandable



PG VR2 (4-key) PG VR2 (6-key) PG VR2 (9-key)

Expandable



Q-Logic 2 EX Joystick Q-Logic 2 EX Joystick with Lights Q-Logic 2 EX Enhanced Display Q-Logic 2 EX Stand Alone Joystick / EX Attendant Joystick Q-Logic 2 Specialty Controls Interface Module

Multi-Switch Systems



ASL 3-Switch Head Array Switch-It Head Array (3, 4 or 5 switch) 3-Switch Stealth™ Ultra Head Array Stealth™ i-Drive® Tri-Array

Proportional Drive Controls



ASL Micro Extremity Control Switch-It Proportional Joystick Stealth™ Mushroom Joystick ASL Rim Control Switch-It Touch Drive

Single Switch



TASH Pillow Switch TASH Trigger Switch Buddy Button Egg Switch MicroLite Switch Switch-It Uniwobble Switch

Sip-N-Puff Options



Therafin Whisper-Lite Sip-N-Puff Hardware ASL Sip-N-Puff Head Array Switch-It Sip-N-Puff with Opti-Stop

Additional Options



ASL Remote Stop Switch Switch-It Remote Stop Switch ASL Remote Attendant Control Joystick Handles

Q-Logic 2 Specialty Controls Interface Module*

The Q-Logic 2 Specialty Controls Interface Module offers an enhanced display at a lower cost and is an economical solution to provide specialty controls. On this module, specialty controls run standard through the Q-Logic 2 Drive Control System. Bluetooth® functionality also operates through Q-Logic 2. The module mounts to the TRU-Balance® 3 Power Positioning Systems accessory bar.

* Either the Q-Logic 2 Enhanced Display or Q-Logic Enhanced Display is still required for Sip-N-Puff and single switch driving.

