



QUANTUM® OUTBACK® 4x4
shown in *Tan* with 10" (25.4cm)
safe seat elevation and 14"
(38.1cm) pneumatic tires



The Quantum® Outback® 4x4 Powerchair is your ultimate adventure companion, ready to tackle almost any terrain with confidence. Designed for superior performance, this rugged powerchair features powerful 4-pole sealed motors and all-wheel drive for unbeatable traction and stability. Its robust construction ensures a smooth and comfortable ride on uneven surfaces, featuring advanced (SRS) suspension and off-road pneumatic 14" knobby tyres. Whether you're navigating gravel paths, bush tracks, or country properties, the Outback® 4x4 powerchair makes it a breeze. Experience unmatched freedom and independence with a powerchair that combines strength, durability, and cutting-edge technology to conquer any challenge!

Q
QUANTUM®
#1 FOR REHAB POWER

quantumrehab.com.au

Specifications

DRIVE WHEELS	Standard: 14" knobby pneumatic (38.1cm)
DRIVETRAIN	4 high torque motors
SUSPENSION	Automotive-grade smooth ride suspension
MAXIMUM SPEED ¹	Up to 10 km/h (5.6 km/h when elevated)
MAXIMUM OBSTACLE CLIMBING ABILITY	10.2 cm
MAXIMUM RATED SLOPE	8° (14.1%)
GROUND CLEARANCE	10.7 cm (at motors)
TURNING RADIUS ²	125 cm
OVERALL LENGTH ²	106.7 cm
BASE WIDTH	73.2 cm
SEATING SIZES	
TRU-Balance® 3 Power Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
TRU-Balance® 3 Power Lift & Tilt	Width: 16-24"(40.6cm - 61cm) Depth: 16-24"(40.6cm - 61cm)
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Power Tilt ^{5,9}	19-20" (48.3 cm - 50.8 cm)
TRU-Balance® 3 Power Lift & Tilt ^{5,9}	20-21" (50.8 cm - 53.4 cm)
BATTERY SIZE	Group 27
BATTERY WEIGHT ⁷	26.6 kg each
RANGE PER CHARGE ³	Up to 25.75 km at 136 kg Up to 27.4 km at 99.8 kg
AVAILABLE ELECTRONICS	120A Q-Logic 3
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-pole motors
WEIGHT CAPACITY ⁴	159 kg
BASE WEIGHT	92 kg

1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
2) Without front riggings
3) Battery capacity may change when adding power positioning systems.
4) Weight capacity is dictated by the selection of seating.
5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 300 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details.
Listed measurements are in the lowest position; does not include elevating system unless specified.
6) Battery weight may vary +/- 3 lbs. based on manufacturer.
7) All Quantum Rehab® units can be ordered with specialty controls.
8) Call for availability
9) Seat-to-floor heights are measured with the unit unloaded.

Standard Features

- 4-pole sealed motors
- Maximum speed of up to 10km/h 5.6 km/h when elevated)
- LED lighting package
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- Unmatched driving performance and ride quality
- Advanced stability for smooth obstacle transitions

Options

- 10" (25.4 cm) of power adjustable safe seat elevation at speeds of up to 5.6 km/h
- Quantum® backup camera
- Glove box
- Backpack holder
- Multi-use accessory bag
- Hydration system
- Accessory bar

Colors

