

NEXT LEVEL REAR-WHEEL DRIVE

The R-TRAK revolutionizes the real-wheel drive power chair experience. Independent suspension on all four wheels combined with a compact base provide a smoother, more manoeuvrable driving experience. Standard features include fender lights, USB charger and drive wheel color accents. The R-TRAK is also available with our 10" safe seat elevation, allowing the client to drive up to 5.6 km/h while elevated.



quantumrehab.com.au

Specifications

DRIVE WHEELS	Standard: 35.5 cm (14") pneumatic Optional: 35.5 cm (14") flat free 35.5 cm (14") knobby
CASTER WHEELS	Front: 22.8 cm (9") solid
ANTI-TIP WHEELS	Rear. 10 cm (4") flat free
MAXIMUM SPEED ¹	Up to 10 km/h (5.6 km/h when elevated)
GROUND CLEARANCE	10 cm (battery tray)
TURNING RADIUS ²	65.2 cm
OVERALL LENGTH ²	88.3 cm
BASE WIDTH	61.4 cm
SEATING SIZES	
TRU-Balance® 3 Static Seat	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Tilt	Width: 12-22" Depth: 12-22"
TRU-Balance® 3 Power Lift & Tilt	Width: 12-22" Depth: 12-22"
SEAT-TO-FLOOR HEIGHTS	
TRU-Balance® 3 Static Seat ⁵	40.6 cm - 20.32 cm (16-18") (not available with Synergy®)
TRU-Balance® 3 Power Tilt ⁵	43.8 cm - 48.8 cm (17.25-19.25")
TRU-Balance® 3 Power Lift & Tilt 5	43.8 cm - 48.8 cm (17.25-19.25")
TRU-Balance® Flex Seating	41.9 cm - 44.4 cm (16.5-17.5")
TRU-Balance® Flex Tilt	41.9 cm (16.25")
BATTERY SIZE ³	Crown 24 (2)
	Group 24 (2)
BATTERY WEIGHT ⁷	Group 24: 23.3 kg each
AVAILABLE ELECTRONICS 8	120A Q-Logic 3
BATTERY CHARGER	8A Off-board
MOTOR PACKAGES	4-pole motors
WEIGHT CAPACITY ⁴	136 kg
BASE WEIGHT	58 kg

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pres-
- 2) Without front riggings
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
 5) All measurements +/- 3 cm with 45.7 cm x 45.7 cm seat, fixed towers, 136 kg. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Battery weight may vary +/- 1.3 kg based on manufacturer.
- 7) All Quantum Rehab® units can be ordered with specialty controls
- 8) Call for availability
- 9) Seat-to-floor heights are measured with the unit unloaded. Due to the R-Trak's suspension system, seat-to-floor heights are approximately 0.63 cm lower when the unit is loaded.

Standard Features

- 4-pole motors
- Maximum speed of up to 10 km/h (5.6 km/h when elevated)
- LED lighting in front caster arms and fenders
- Superior, automotive-grade SRS (Smooth Ride Suspension)
- · Unmatched driving performance and ride quality
- · Advanced stability for smooth obstacle transitions

Options

- · 25.4 cm (10") of power adjustable safe seat elevation at speeds of up to 5.6 km/h
- · Quantum® backup camera
- Trilogy vent tray 8
- Compact vent tray ⁸

Colours





Midnight Blue



Electric Blue



Back in Black

Night Train



Lemon Crush







